

Product Datasheet

Anti-Human CD62L Monoclonal Antibody, Azide Free Clone B-S13 852.780

Product Name Anti-Human CD62L Monoclonal Antibody, Azide Free Clone B-S13

Catalog Number 852.780

Description Mouse monoclonal

Species/Host Mouse

Conjugation Azide free

Tested Applications ELISA, FC

Immunogen CD62L transfected CHO cell line

Brand Diaclone

Form/Appearance Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free

Storage Stable at +2-8°C for 12 months. For longer storage freeze aliquots.

Isotype Mouse IgG1 Kappa light chain

Clonality Monoclonal

Specificity Recognises the Leucocyte Endothelial Cell Adhesion Molecule-1 (LECAM-1), L-selectin, a 74-95 kDa protein. Cross reacts with mouse CD62L protein in cytometry.

Hybridoma Myeloma X63/AG.8653 x Balb/c spleen cells