

Product Datasheet

Anti-Human CD62L Monoclonal Antibody, Unconjugated Clone B-S13 852.783.020

Product Name	Anti-Human CD62L Monoclonal Antibody, Unconjugated Clone B-S13
Catalog Number	852.783.020
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Purified
Tested Applications	FC
Immunogen	CD62L transfected CHO cell line
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline with 1% BSA and 0.09% sodium azide
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Note	Use 10 µl to label 10 ⁶ cells or 100 µl of whole blood
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises the Leucocyte Endothelial Cell Adhesion Molecule-1 (LECAM-1), L-selectin, a 74-95 kDa protein; Cross reacts with mouse CD62L protein in cytometry.
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells