

## Product Datasheet

### Anti-Human Galectin-10 Monoclonal Antibody, Unconjugated Clone B-F42 852.963.020

Product Name	Anti-Human Galectin-10 Monoclonal Antibody, Unconjugated Clone B-F42
Catalog Number	852.963.020
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Purified
Tested Applications	FC, IHC
Immunogen	Recombinant human Galectin-10
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline with 1% BSA and 0.09% sodium azide
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Note	Galectin-10 is an intracellular antigen. Fixation with paraformaldehyde 4% and permeabilisation with saponin 0.1% are required before staining. Use 10 µl to label 10 <sup>6</sup> cells
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises both natural and recombinant human Galectin-10
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells