

Product Datasheet

Anti-Human CD261 Monoclonal Antibody, Azide Free Clone B-N28 852.980

Product Name Anti-Human CD261 Monoclonal Antibody, Azide Free Clone B-N28

Catalog Number 852.980

Description Mouse monoclonal

Species/Host Mouse

Conjugation Azide free

Tested Applications IP, WB

Immunogen Recombinant human TRAIL R1/Fc chimera

Brand Diaclone

Form/Appearance Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free

Storage Stable at +2-8°C for 12 months

Isotype Mouse IgG1 Kappa light chain

Clonality Monoclonal

Specificity Recognises a 43-46 kDa protein

Hybridoma Myeloma X63/AG.8653 x Balb/c spleen cells