

Product Datasheet

Anti-Human p53 Monoclonal Antibody, Azide Free Clone B-P3 853.040

Product Name	Anti-Human p53 Monoclonal Antibody, Azide Free Clone B-P3
Catalog Number	853.040
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA, IP, WB, IHC, FC
Immunogen	Recombinant p53 protein
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG2a Kappa light chain
Clonality	Monoclonal
Specificity	Recognises wild type and mutant p53, cross-react with p53 in other primates. The epitope is within residues 18-30 of human p53
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells