

Product Datasheet

Anti-Human CD253 Monoclonal Antibody, Azide Free Clone B-T24 853.090

Product Name Anti-Human CD253 Monoclonal Antibody, Azide Free Clone B-T24

Catalog Number 853.090

Description Mouse monoclonal

Species/Host Mouse

Conjugation Azide free

Tested Applications ELISA, Functional assay, FC

Immunogen Recombinant human TRAIL

Brand Diaclone

Form/Appearance Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free

Storage Stable at +2-8°C for 12 months. For longer storage freeze aliquots.

Isotype Mouse IgG1 Kappa light chain

Clonality Monoclonal

Specificity Recognises the TRAIL, APO2L antigen, a 32 kDa protein

Biological Activity Inhibits TRAIL induced apoptosis

Hybridoma Myeloma X63/AG.8653 x Balb/c spleen cells