

Product Datasheet

Alpha-2 Antiplasmin, Human Plasma 107-50

Product Name	Alpha-2 Antiplasmin, Human Plasma
Catalog Number	107-50
Alternative Names	A2AP, alpha 2-antiplasmin, Plasmin Inhibitor
Brand	Lee Biosolutions
Form/Appearance	Lyophilized from 20 mM Bis-Tris, pH 6.4, with 200 mM NaCl.
Concentration	Reported in mg/vial (A_{280} , $E^{0.1\%}=0.7$)
Purity	$\geq 95\%$ (SDS-PAGE)
Source	Human Plasma
Molecular Weight	63-70 kDa
Shelf Life	2 years
Specific Activity	≥ 5 IU/mg protein