

## Product Datasheet

### Apolipoprotein CIII (ApoCIII), Human - Ultra Pure 125-21

**Product Name** Apolipoprotein CIII (ApoCIII), Human - Ultra Pure

**Catalog Number** 125-21

**Alternative Names** ApoCIII, ApoC3

**Brand** Lee Biosolutions

**Form/Appearance** Liquid

**Concentration** Typically 1.0 mg/mL (Lowry)

**Storage** -20 °C

**Purity** ≥ 98% (SDS-PAGE)

**Source** Human Plasma

**Shelf Life** 2 years