

Product Datasheet

C-Reactive Protein (CRP), Human Pleural Fluid 140-11

Product Name C-Reactive Protein (CRP), Human Pleural Fluid

Catalog Number 140-11

Brand Lee Biosolutions

Form/Appearance Liquid Clear solution

Concentration 1.0 - 5.0 mg/mL (A_{280nm} , $E^{0.1\%} = 1.8$)

Storage 2-8° C

Purity ≥ 99% (SDS-PAGE)

CAS Number 9007-41-4

Source Human Pleural Fluid

Molecular Weight 118000

Shelf Life 3 years