

Product Datasheet

Anti-Human CD1a Monoclonal Antibody, Unconjugated Clone HI149 853.153.020

Product Name	Anti-Human CD1a Monoclonal Antibody, Unconjugated Clone HI149
Catalog Number	853.153.020
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Purified
Tested Applications	IHC, FC
Immunogen	Human thymocytes
Brand	Diaclone
Form/Appearance	Concentration of 1mg/1ml in phosphate buffered saline containing 15mM sodium azide
Storage	Store at +2-8°C. DO NOT FREEZE.
Note	We recommend that each laboratory determine an optimum working titre for use in its particular application
Isotype	IgG1
Clonality	Monoclonal
Specificity	Reacts with CD1a (T6), a 49 KDa polypeptide associated with b2-microglobulin expressed on cortical thymocytes (strongly), Langerhans cells and dendritic cells