

Product Datasheet

Recombinant Human C-Reactive Protein (rCRP) 140-11R

Product Name	Recombinant Human C-Reactive Protein (rCRP)
Catalog Number	140-11R
Alternative Names	r-CRP, r-hCRP
Brand	Lee Biosolutions
Form/Appearance	Liquid
Concentration	0.5 - 2.0 mg/mL ($E^{0.1\%}_{280nm} = 1.8$)
Storage	2-8° C
Purity	Highly Purified
CAS Number	99401-15-7
Source	Recombinant
Buffer	20 mM Tris HCl, 140 mM NaCl, 2 mM CaCl ₂ , with 0.05% sodium azide, pH 7.5.
Molecular Weight	23047
Shelf Life	4 years
Preservatives	0.05% sodium azide