

Product Datasheet

C-Reactive Protein (CRP), Canine (Dog) 143-10

Product Name C-Reactive Protein (CRP), Canine (Dog)

Catalog Number 143-10

Alternative Names Dog CRP, Canine CRP

Brand Lee Biosolutions

Form/Appearance Liquid

Concentration ≥ 0.5 mg/mL ($E^{1\text{ mg/mL}}_{280\text{nm}} = 1.6$)

Storage -60 °C or below

Purity $\geq 95\%$ (SDS-PAGE)

Source Canine Serum

Molecular Weight $\sim 20,000$

Shelf Life 2 years