

## Product Datasheet

### Carcinoembryonic Antigen (CEA), Human - High Purity 151-10

<b>Product Name</b>	Carcinoembryonic Antigen (CEA), Human - High Purity
<b>Catalog Number</b>	151-10
<b>Alternative Names</b>	CEA, Carcinoembryonic antigen-related cell adhesion molecule 5, Meconium antigen 100, CD_antigen=CD66e
<b>Brand</b>	Lee Biosolutions
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	Report mg/mL (A <sub>280</sub> and Roche EIA)
<b>Storage</b>	-20 °C
<b>Purity</b>	≥ 95% (SDS-PAGE)
<b>Source</b>	Human Metastatic Liver Carcinoma
<b>Molecular Weight</b>	~180,000
<b>Shelf Life</b>	4 years