

Product Datasheet

Carcinoembryonic Antigen (CEA), Human - High Purity 151-10

Product Name	Carcinoembryonic Antigen (CEA), Human - High Purity
Catalog Number	151-10
Alternative Names	CEA, Carcinoembryonic antigen-related cell adhesion molecule 5, Meconium antigen 100, CD_antigen=CD66e
Brand	Lee Biosolutions
Form/Appearance	Liquid
Concentration	Report mg/mL (A ₂₈₀ and Roche EIA)
Storage	-20 °C
Purity	≥ 95% (SDS-PAGE)
Source	Human Metastatic Liver Carcinoma
Molecular Weight	~180,000
Shelf Life	4 years