

Product Datasheet

Creatine Kinase MB (rCK-MB), Recombinant Human 190-24-R1

Product Name	Creatine Kinase MB (rCK-MB), Recombinant Human
Catalog Number	190-24-R1
Tested Applications	ELISA, WB
Brand	Lee Biosolutions
Form/Appearance	Liquid Clear, colorless solution
Concentration	≥ 1.0 mg/mL (Coomassie)
Storage	-70 °C or below
Purity	≥ 95% (SDS-PAGE)
Source	Recombinant
Shelf Life	2 years
Biological Activity	≥ 400 U/mL @ 37°C
Unit Definition	One unit will catalyze the transphosphorylation of one micromole of phosphate from creatine phosphate to ADP per minute at 37°C.