

Product Datasheet

Recombinant Human Creatine Kinase MB (rCK-MB) 190-24-R2

Product Name	Recombinant Human Creatine Kinase MB (rCK-MB)
Catalog Number	190-24-R2
Alternative Names	r-hCKMB II, r-hCK-MB, r-hCK-2
Brand	Lee Biosolutions
Form/Appearance	Liquid
Concentration	≥ 10 mg/mL (Coomassie)
Storage	-70 °C or below
Purity	Highly Purified
CAS Number	9001-15-4
Source	Recombinant
Molecular Weight	84000
Shelf Life	5 years
Biological Activity	≥ 10,000 U/mL
Unit Definition	One unit will catalyze the transphosphorylation of one micromole of phosphate from creatine phosphate to ADP per minute at 37°C.