

## Product Datasheet

### Recombinant Human Creatine Kinase MB (rCK-MB) 190-24-R2

<b>Product Name</b>	Recombinant Human Creatine Kinase MB (rCK-MB)
<b>Catalog Number</b>	190-24-R2
<b>Alternative Names</b>	r-hCKMB II, r-hCK-MB, r-hCK-2
<b>Brand</b>	Lee Biosolutions
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	≥ 10 mg/mL (Coomassie)
<b>Storage</b>	-70 °C or below
<b>Purity</b>	Highly Purified
<b>CAS Number</b>	9001-15-4
<b>Source</b>	Recombinant
<b>Molecular Weight</b>	84000
<b>Shelf Life</b>	5 years
<b>Biological Activity</b>	≥ 10,000 U/mL
<b>Unit Definition</b>	One unit will catalyze the transphosphorylation of one micromole of phosphate from creatine phosphate to ADP per minute at 37°C.