

Product Datasheet

Creatine Kinase MM (CKMM) Goat anti-Human Polyclonal Antibody 190-46

Product Name Creatine Kinase MM (CKMM) Goat anti-Human Polyclonal Antibody

Catalog Number 190-46

Species/Host Goat

Tested Applications ELISA, Blotting

Immunogen High purity human Creatine Kinase MM (CKMM)

Brand Lee Biosolutions

Form/Appearance Liquid

Concentration ≥ 1.0 mg/mL ($E^{0.1\%}_{280nm} = 1.4$)

Storage -20 °C

Purity Affinity Purified

Source Polyclonal Antibody

Shelf Life 2 years