

---

## Product Datasheet

### D-Dimer from human plasma 200-12

|                          |   |
|--------------------------|---|
| <b>Product Name</b>      | D-Dimer from human plasma                     |
| <b>Catalog Number</b>    | 200-12  |
| <b>Alternative Names</b> | Fragment D-dimer, Fibrin degradation fragment |
| <b>Brand</b>             | Lee Biosolutions                              |
| <b>Form/Appearance</b>   | Liquid Clear, colorless solution              |
| <b>Storage</b>           | -20 °C  |
| <b>Purity</b>            | ≥ 80% (SDS-PAGE)                              |
| <b>Source</b>            | Human Plasma                                  |
| <b>Molecular Weight</b>  | ~200,000                                      |
| <b>Shelf Life</b>        | 2 years                                       |