

Product Datasheet

D-Dimer from human plasma 200-13

Product Name	D-Dimer from human plasma
Catalog Number	200-13
Alternative Names	Fragment D-dimer, Fibrin degradation fragment
Brand	Lee Biosolutions
Form/Appearance	Liquid Clear, colorless solution
Concentration	1 - 5 mg/mL ($E_{0.1\%}^{0.1\%280nm} = 1.78$)
Storage	-20 °C
Purity	≥ 95% (SDS-PAGE)
Source	Human Plasma
Molecular Weight	~200,000
Shelf Life	2 years