

Product Datasheet

D-Dimer, anti-Human Monoclonal Antibody 200-14

Product Name D-Dimer, anti-Human Monoclonal Antibody

Catalog Number 200-14

Species/Host Mouse

Immunogen Purified, native human D-Dimer

Alternative Names anti-D-dimer

Brand Lee Biosolutions

Form/Appearance Liquid

Concentration Reported in mg/mL ($E^{0.1\%}_{280nm} = 1.4$)

Storage -20 °C

Isotype IgG

Purity Highly Purified

Source Monoclonal Antibody

Shelf Life 2 years