

## Product Datasheet

### Ferritin, Human Liver 270-40

<b>Product Name</b>	Ferritin, Human Liver
<b>Catalog Number</b>	270-40
<b>Alternative Names</b>	FTL
<b>Brand</b>	Lee Biosolutions
<b>Form/Appearance</b>	Solution in tris buffered sodium chloride, pH 7.5, with bromo-nitro-dioxane and methylisothiazolone as preservatives. Clear, reddish-brown solution
<b>Concentration</b>	3 - 6 mg/mL (Lowry)
<b>Storage</b>	2-8° C
<b>Purity</b>	≥ 95% (SDS-PAGE)
<b>CAS Number</b>	9007-73-2
<b>Source</b>	Human Liver
<b>Molecular Weight</b>	440000
<b>Shelf Life</b>	3 years