

## Product Datasheet

### Ferritin, Human Placenta 270-60

<b>Product Name</b>	Ferritin, Human Placenta
<b>Catalog Number</b>	270-60
<b>Brand</b>	Lee Biosolutions
<b>Form/Appearance</b>	Liquid Clear, reddish-brown solution
<b>Concentration</b>	2.0 - 7.0 mg/mL (Lowry)
<b>Storage</b>	2-8° C
<b>Purity</b>	≥ 96% (SDS-PAGE)
<b>CAS Number</b>	9007-73-2
<b>Source</b>	Human Placenta
<b>Molecular Weight</b>	440000
<b>Shelf Life</b>	2 years