

Product Datasheet

Galectin-3 antigen 280-15

Product Name	Galectin-3 antigen
Catalog Number	280-15
Alternative Names	GAL-3, LGALS3
Brand	Lee Biosolutions
Form/Appearance	Liquid
Concentration	≥ 0.5 mg/mL (Immunoassay)
Storage	-20 °C
Purity	≥ 90% (SDS-PAGE)
Source	Recombinant
Molecular Weight	~30,000
Shelf Life	1 year
Biological Activity	Immunologically active