

Product Datasheet

Kallikrein, Porcine Pancreas 314-10

Product Name	Kallikrein, Porcine Pancreas
Catalog Number	314-10
Alternative Names	Kininogenase, Kininogenin
Brand	Lee Biosolutions
Form/Appearance	Lyophilized White to brown powder
Purity	Purified
CAS Number	2593801
Source	Porcine Pancreas
Shelf Life	2 years
Biological Activity	≥ 50 U/mg
Specific Activity	≥ 300 U/mg protein
Unit Definition	One unit of Kallikrein will hydrolyze 1.0 μ mole of N α -benzoyl-L-arginine ethyl ester (BASE) to N α -benzoyl-L-arginine and ethanol per minute at pH 8.0 at 25°C.