

## Product Datasheet

### Haptoglobin (Hp) - Mixed Genotype from human plasma 330-11

**Product Name** Haptoglobin (Hp) - Mixed Genotype from human plasma

**Catalog Number** 330-11

**Alternative Names** HAP, HP, Hp 1-1, 2-1, and 2-2

**Brand** Lee Biosolutions

**Form/Appearance** Lyophilized White to near white

**Concentration** Reported in mg/vial ( $A_{278} E^{0.1\%} = 1.20$ )

**Purity**  $\geq 95\%$  (SDS-PAGE)

**CAS Number** 9087-69-8

**Source** Human Plasma

**Molecular Weight** 86,000, 200,000 and 400,000

**Shelf Life** 5 years