

Product Datasheet

Haptoglobin (Hp) - Mixed Genotype from human plasma 330-11

Product Name Haptoglobin (Hp) - Mixed Genotype from human plasma

Catalog Number 330-11

Alternative Names HAP, HP, Hp 1-1, 2-1, and 2-2

Brand Lee Biosolutions

Form/Appearance Lyophilized White to near white

Concentration Reported in mg/vial ($A_{278} E^{0.1\%} = 1.20$)

Purity $\geq 95\%$ (SDS-PAGE)

CAS Number 9087-69-8

Source Human Plasma

Molecular Weight 86,000, 200,000 and 400,000

Shelf Life 5 years