

Product Datasheet

Haptoglobin, Mouse Serum 334-11

Product Name Haptoglobin, Mouse Serum

Catalog Number 334-11

Alternative Names Hp

Brand Lee Biosolutions

Form/Appearance Liquid Frozen solution

Concentration ≥ 0.5 mg/mL (Coomassie)

Storage -20 °C

Purity $\geq 90\%$ (SDS-PAGE)

CAS Number 9087-69-8

Source Mouse Serum

Molecular Weight alpha Hp: $\sim 10,000$
beta Hp: $\sim 40,000$

Shelf Life 2 years