

## Product Datasheet

### Immunoglobulin A (IgA), Human Plasma 340-10

<b>Product Name</b>	Immunoglobulin A (IgA), Human Plasma
<b>Catalog Number</b>	340-10
<b>Alternative Names</b>	IgA
<b>Brand</b>	Lee Biosolutions
<b>Form/Appearance</b>	0.2 micron filtered solution in 10 mM sodium phosphate, pH 7.4, containing 150 mM sodium chloride and 0.09% sodium azide. Clear, colorless to slight yellow solution
<b>Concentration</b>	3 - 6 mg/mL ( $A_{280}$ , $E^{0.1\%} = 1.32$ )
<b>Storage</b>	2-8° C
<b>Purity</b>	≥ 98% (SDS-PAGE)
<b>Source</b>	Human Plasma
<b>Specificity</b>	Mucosal Immunity
<b>Molecular Weight</b>	160000
<b>Shelf Life</b>	2 years