

Product Datasheet

Anti-Human CD36 Monoclonal Antibody, Unconjugated Clone TR9 853.213.020

Product Name Anti-Human CD36 Monoclonal Antibody, Unconjugated Clone TR9

Catalog Number 853.213.020

Description Mouse monoclonal

Species/Host Mouse

Conjugation Purified

Tested Applications FC

Immunogen Platelets

Brand Diaclone

Form/Appearance 20 mM Tris/HCl with 15 mM sodium azide, approx. pH 7.4

Storage Store at +2-8°C. DO NOT FREEZE.

Note We recommend that each laboratory determine an optimum working titre for use in its particular application

Isotype IgG1

Clonality Monoclonal

Source Platelets

Specificity Reacts with CD36 (GPIIb), a 85 kDa integral membrane glycoprotein expressed on platelets, macrophages, endothelial cells, early erythroid cells and megakaryocytes