

Product Datasheet

Immunoglobulin A (IgA), Human Plasma 340-12

Product Name Immunoglobulin A (IgA), Human Plasma

Catalog Number 340-12

Alternative Names IgA

Brand Lee Biosolutions

Form/Appearance Lyophilized White to off-white powder

Concentration Reported in mg/vial (A_{280} , $E_{1\text{mg/mL}} = 1.34$)

Storage -20 °C

Purity $\geq 96\%$ (SDS-PAGE)

Source Human Plasma

Shelf Life 5 years