

Product Datasheet

Immunoglobulin A (IgA), Human Plasma 340-13

Product Name	Immunoglobulin A (IgA), Human Plasma
Catalog Number	340-13
Alternative Names	IgA, slgA
Brand	Lee Biosolutions
Form/Appearance	Liquid Clear, colorless to slight yellow solution
Concentration	3 - 6 mg/mL (A ₂₈₀ , E ^{0.1%} = 1.32)
Storage	2-8° C
Purity	≥ 95% (SDS-PAGE)
Source	Human Plasma
Specificity	Mucosal Immunity
Molecular Weight	160000
Shelf Life	2 years