

Product Datasheet

Immunoglobulin G Fab (IgG Fab) Antigen from human plasma 340-27

Product Name Immunoglobulin G Fab (IgG Fab) Antigen from human plasma

Catalog Number 340-27

Alternative Names IgG Fab, Fragment antigen-binding

Brand Lee Biosolutions

Form/Appearance Liquid Clear, colorless liquid

Concentration ≥ 5.0 mg/mL ($E^{1\text{ mg/mL}}_{280\text{nm}} = 1.35$)

Storage -20 °C

Purity $\geq 95\%$ (SDS-PAGE)

Source Human Plasma

Molecular Weight $\sim 50,000$

Shelf Life 3 years