

## Product Datasheet

### Immunoglobulin M (IgM), Human Plasma 340-30

<b>Product Name</b>	Immunoglobulin M (IgM), Human Plasma
<b>Catalog Number</b>	340-30
<b>Alternative Names</b>	IgM
<b>Brand</b>	Lee Biosolutions
<b>Form/Appearance</b>	Lyophilized from 50 mM tris, 200 mM sodium chloride, and 0.05% sodium azide, pH 8.0. White to off-white powder
<b>Concentration</b>	Reported in mg/vial ( $E^{0.1\%}_{280nm} = 1.18$ )
<b>Storage</b>	2-8° C
<b>Purity</b>	≥ 96% (SDS-PAGE)
<b>Source</b>	Human Plasma
<b>Molecular Weight</b>	950000
<b>Shelf Life</b>	5 years