

Product Datasheet

Immunoglobulin M (IgM), Human Plasma 340-31

Product Name	Immunoglobulin M (IgM), Human Plasma
Catalog Number	340-31
Alternative Names	IgM
Brand	Lee Biosolutions
Form/Appearance	Frozen solution in 50 mM tris, 200 mM sodium chloride, and 0.05% sodium azide, pH 8.0.
Concentration	1.0 - 3.0 mg/mL ($E^{0.1\%}_{280nm} = 1.18$)
Storage	-20 °C
Purity	≥ 95% (SDS-PAGE)
Source	Human Plasma
Molecular Weight	950000
Shelf Life	5 years