

## Product Datasheet

### Immunoglobulin E (IgE) - Human Myeloma Plasma 340-40

<b>Product Name</b>	Immunoglobulin E (IgE) - Human Myeloma Plasma
<b>Catalog Number</b>	340-40
<b>Alternative Names</b>	IgE
<b>Brand</b>	Lee Biosolutions
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	-20 °C
<b>Source</b>	Human Myeloma Plasma
<b>Shelf Life</b>	Due to differing stabilities of various biomarkers in this complex biological material it is not possible to assign a single expiration date that is relevant to all biomarkers. The end user should monitor the biomarkers of interest to verify suitability over time.
<b>Biological Activity</b>	≥ 1,500,000 U/mL (ELISA)