

## Product Datasheet

### Interleukin-6 (IL-6) Antigen 341-36

<b>Product Name</b>	Interleukin-6 (IL-6) Antigen
<b>Catalog Number</b>	341-36
<b>Alternative Names</b>	IL6, BSF2, HGF, HSF, IFNB2, IL-6, BSF-2, CDF, IFN-beta-2, interleukin 6
<b>Brand</b>	Lee Biosolutions
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	≥ 0.5 mg/mL (Immunoassay)
<b>Storage</b>	-20 °C
<b>Purity</b>	≥ 90% (SDS-PAGE)
<b>Source</b>	Recombinant
<b>Molecular Weight</b>	~20,813 (excluding tag)
<b>Shelf Life</b>	1 year
<b>Biological Activity</b>	Immunologically active