

Product Datasheet

Recombinant NGAL Antigen 342-49R

| | |
|--------------------------|--|
| Product Name | Recombinant NGAL Antigen |
| Catalog Number | 342-49R |
| Alternative Names | NGAL, LCN2, Lipocalin-2, oncogene 24p3 |
| Brand | Lee Biosolutions |
| Form/Appearance | Lyophilized powder from 0.02 M ammonium bicarbonate. |
| Concentration | Reported in mg/vial (BCA) |
| Storage | -20 °C |
| Purity | ≥ 96% (SDS-PAGE) |
| Source | Recombinant |
| Molecular Weight | ~21,000 |
| Shelf Life | 2 years |