

## Product Datasheet

### Low Density Lipoprotein (LDL) - Low Triglycerides, Human 360-13

<b>Product Name</b>	Low Density Lipoprotein (LDL) - Low Triglycerides, Human
<b>Catalog Number</b>	360-13
<b>Alternative Names</b>	LDL, LDL-C, Beta-Lipoprotein, Concentrate, Fraction
<b>Brand</b>	Lee Biosolutions
<b>Form/Appearance</b>	Liquid Yellow to orange liquid
<b>Storage</b>	-60 °C or below
<b>Purity</b>	Plasma Fraction
<b>Source</b>	Human Plasma
<b>Molecular Weight</b>	2.75 million (ref. Tetz second Edition)
<b>Shelf Life</b>	2 years
<b>Content</b>	4,000 - 6,000 mg/dL (Enzymatic) <b>Inquire for custom concentrations and cholesterol ratios</b>
<b>Triglyceride Ratio</b>	≤ 0.80 (Triglycerides/Total Cholesterol)