

Product Datasheet

Low Density Lipoprotein (LDL) - Preservative Added 360-16

| | |
|---------------------------|--|
| Product Name | Low Density Lipoprotein (LDL) - Preservative Added |
| Catalog Number | 360-16 |
| Alternative Names | LDL, LDL-C, Beta-Lipoprotein, Concentrate, Fraction |
| Brand | Lee Biosolutions |
| Form/Appearance | Liquid Yellow to orange liquid |
| Storage | 2-8° C |
| Purity | Plasma Fraction |
| Source | Human Plasma |
| Molecular Weight | 2.75 million (ref. Tetz second Edition) |
| Shelf Life | 1 year |
| Content | ≥ 3,000 mg/dL (Enzymatic) Inquire for custom concentrations and cholesterol ratios |
| Triglyceride Ratio | ≤ 0.9 (Triglycerides/Total Cholesterol) |