

## Product Datasheet

### anti-Human Lipoprotein(a) [Lp(a)] Polyclonal Antibody 400-45

<b>Product Name</b>	anti-Human Lipoprotein(a) [Lp(a)] Polyclonal Antibody
<b>Catalog Number</b>	400-45
<b>Species/Host</b>	Sheep
<b>Immunogen</b>	High purity human Lp(a)
<b>Alternative Names</b>	Lp(a), LPA, AK38, APOA, LP
<b>Brand</b>	Lee Biosolutions
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	0.5 - 2.0 mg/mL (A <sub>280nm</sub> )
<b>Storage</b>	-20 °C
<b>Purity</b>	Affinity Purified
<b>Source</b>	Polyclonal Antibody
<b>Shelf Life</b>	1 year