

Product Datasheet

Lipoprotein-associated Phospholipase A2 (Lp-PLA2) antigen 400-60

Product Name Lipoprotein-associated Phospholipase A2 (Lp-PLA2) antigen

Catalog Number 400-60

Alternative Names Lp-PLA2, platelet-activating factor acetylhydrolase, PAF acetylhydrolase, PLA2G7

Brand Lee Biosolutions

Form/Appearance Liquid

Concentration ≥ 0.5 mg/mL

Storage -20 °C

Purity $\geq 90\%$ (SDS-PAGE)

Source Recombinant

Molecular Weight 45000

Shelf Life 18 months

Biological Activity Immunologically active