

## Product Datasheet

### Kappa Free Light Chain antigen from human urine 410-12

**Product Name** Kappa Free Light Chain antigen from human urine

**Catalog Number** 410-12

**Alternative Names** Bence Jones Protein-Kappa, FLC-Kappa, KLC, KFLC

**Brand** Lee Biosolutions

**Form/Appearance** Liquid Colorless to pale yellow

**Concentration**  $\geq 1.0$  mg/mL (Nephelometry Method)

**Storage** 2-8° C

**Purity**  $\geq 95\%$  (SDS-PAGE)

**Source** Human Urine

**Molecular Weight** ~23,000

**Shelf Life** 18 months

**Content** Lambda Free Light Chain:  $\leq 1\%$   
IgG:  $\leq 1\%$   
IgM:  $\leq 1\%$   
IgA:  $\leq 1\%$   
Cystatin C:  $\leq 1\%$   
Microalbumin:  $\leq 1\%$