

## Product Datasheet

### Kininogen - High Molecular Weight, Human Plasma 411-15

<b>Product Name</b>	Kininogen - High Molecular Weight, Human Plasma
<b>Catalog Number</b>	411-15
<b>Alternative Names</b>	Williams-Fitzgerald-Flaujeac factor, HMWK-kallikrein factor, HK
<b>Brand</b>	Lee Biosolutions
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	$\geq 1$ mg/mL ( $A_{280nm}$ , $E = 0.70$ )
<b>Purity</b>	$\geq 95\%$ (SDS-PAGE)
<b>Source</b>	Human Plasma
<b>Molecular Weight</b>	120000
<b>Shelf Life</b>	2 years