

Product Datasheet

Kininogen - High Molecular Weight, Human Plasma 411-15

Product Name	Kininogen - High Molecular Weight, Human Plasma
Catalog Number	411-15
Alternative Names	Williams-Fitzgerald-Flaujeac factor, HMWK-kallikrein factor, HK
Brand	Lee Biosolutions
Form/Appearance	Liquid
Concentration	≥ 1 mg/mL (A_{280nm} , $E = 0.70$)
Purity	$\geq 95\%$ (SDS-PAGE)
Source	Human Plasma
Molecular Weight	120000
Shelf Life	2 years