

Product Datasheet

Myoglobin (Mb), Human Heart 431-11

Product Name	Myoglobin (Mb), Human Heart
Catalog Number	431-11
Brand	Lee Biosolutions
Form/Appearance	Solution in phosphate buffered saline, pH 7.5 with 0.09% sodium azide. Clear, reddish-brown solution
Concentration	2 - 4 mg/mL (Lowry)
Storage	2-8° C
Purity	≥ 99% (SDS-PAGE)
CAS Number	11080-17-4
Source	Human Heart
Molecular Weight	17700
Shelf Life	3 years