

## Product Datasheet

### Prealbumin/Transthyretin, Human Plasma 496-11

<b>Product Name</b>	Prealbumin/Transthyretin, Human Plasma
<b>Catalog Number</b>	496-11
<b>Alternative Names</b>	Transthyretin, Thyroxine Binding Prealbumin, ATTR, TBPA, TTR
<b>Brand</b>	Lee Biosolutions
<b>Form/Appearance</b>	Lyophilized White to off-white lyophilized solid
<b>Concentration</b>	Reported in mg protein/vial ( $A_{280}$ , $E^{0.1\%}=1.4$ )
<b>Storage</b>	2-8° C
<b>Purity</b>	≥ 96% (SDS-PAGE)
<b>CAS Number</b>	87090-18-4
<b>Source</b>	Human Plasma
<b>Molecular Weight</b>	55000
<b>Shelf Life</b>	5 years