

Product Datasheet

Anti-Human CD105 Monoclonal Antibody, Unconjugated Clone MEM-226 853.263.020

| | |
|----------------------------|---|
| Product Name | Anti-Human CD105 Monoclonal Antibody, Unconjugated Clone MEM-226 |
| Catalog Number | 853.263.020 |
| Description | Mouse monoclonal |
| Species/Host | Mouse |
| Conjugation | Purified |
| Tested Applications | IP, WB, FC |
| Immunogen | Recombinant Vaccinia virus |
| Brand | Diaclone |
| Form/Appearance | Phosphate buffered saline (PBS) with 15 mM sodium azide, approx. pH 7.4 |
| Storage | Store at +2-8°C. DO NOT FREEZE. |
| Note | We recommend that each laboratory determine an optimum working titre for use in its particular application |
| Isotype | IgG2a |
| Clonality | Monoclonal |
| Specificity | Reacts with CD105 (Endoglin), a 180 kDa type I integral membrane homodimer glycoprotein expressed on vascular endothelial cells |