

Product Datasheet

Procalcitonin (PCT) antigen 496-70

| | |
|----------------------------|----------------------------------|
| Product Name | Procalcitonin (PCT) antigen |
| Catalog Number | 496-70 |
| Alternative Names | PCT, Calcitonin precursor |
| Brand | Lee Biosolutions |
| Form/Appearance | Liquid |
| Concentration | ≥ 0.5 mg/mL |
| Storage | -20 °C |
| Purity | ≥ 90% (SDS-PAGE) |
| Source | Recombinant |
| Molecular Weight | 12,796 (sequence, excluding tag) |
| Shelf Life | 18 months |
| Biological Activity | Immunologically active |