

Product Datasheet

Prostate Specific Antigen (PSA), Human Semen 497-22

Product Name	Prostate Specific Antigen (PSA), Human Semen
Catalog Number	497-22
Alternative Names	Gamma-seminoprotein, Seminol, Kallikrein-3, P-30 antigen, Semenogelase
Brand	Lee Biosolutions
Form/Appearance	Lyophilized White to off-white powder
Concentration	Report mg protein/vial (A_{280} , $E^{0.1\%}=1.7$)
Storage	2-8° C
Purity	≥ 96% (SDS-PAGE)
Source	Human Semen
Molecular Weight	28430
Shelf Life	5 years